



Cable Management

UA1533836



To be mounted sideways at the multi function rack or between both multi function racks; for hiding excess cable between medical device to electrical socket

Weight: 0,4 kg:

Mounting type: Round tube Ø 38 mm

COMPATIBLE WITH

CEILING UNIT M20 W/VALIA X-V - CEILING UNIT M20 W/MEDILIFT-H - CEILING UNIT M20 W/MEDILIFT-V -
CEILING UNIT M20 W/MEDIBOOM-H - CEILING UNIT M20 W/MEDIBOOM-V - CEILING UNIT M20 W/VALIA X-H

