



Crutch Stirrup Pad

BF48



Crutch Stirrup Pad

Dimensions : 39,4 x 24,1 x 1,3 cm (15,5" x 9,5" x 0,5")

The OASIS gel positioners combine lightweight foam with silicone gel and are shaped for specific patient positioning needs

Not CE marked product

COMPATIBLE WITH

CMAX™ X-RAY - 5085 SRT - CMAX™ 3 - CMAX™ 3-DRIVE - KERNA™ - 4095 - 5095 - 7080

