




Compensated Single Articulation Headrest (EU rails)



This Head Rest Extension for the CMAX X-Ray is installed using (2) Head Rest Adaptors (BF790) on either end of the table.

With the ability to easily tilt the patient's head, this headrest creates convenient positioning for easier intubation and facilitates support of the patient's head, eliminating the need for additional positioning devices. The additional 12" of tabletop length increases the available imaging area for the image-guided surgical procedure. Mattress pad not included.

Dimensions: 13" (L) x 19.7" (W). Tilt: 90° / -90°.

Weight: 13.2 lbs. (6 kg)

To be used with EU rails

Replacement Pads: 2" TLT Mattress Pad: BF581 / 3"
TLT Mattress Pad: BF231

COMPATIBLE WITH
CMAX™ X-RAY

