



Joystick Controller for CMAX™ X-Ray (EU Side Rails)

TAB800JS



Through palm contact of the orange joystick, this table operating control generates precise movements and allows the surgeon or OR caregiver to easily maneuver the carbon fiber tabletop without any associated vibration. Can be used in conjunction with the hand control.

The Bluetooth-enabled option eliminates the hanging cord promoting an artifact-free imaging environment.

Dimensions: 5 7/8" (L) x 2 3/4" (W) x 5 5/8" (H). Weight: 2.8 lbs.

To be used with EU rails only

- Bluetooth or wired capability
- Floating control knob
- Float Lock/Unlock
- Height Up/Down
- Auto Level
- Indicative LED

COMPATIBLE WITH
CMAX™ X-RAY

